


0248 TRIMMING INSTRUCTIONS

***** NOTE *****

SOME INCUBATOR TRAYS ARE TOO SHORT TO ACCEPT TWO RACKS.

WHEN THIS HAPPENS, USE A KNIFE, HACK SAW OR BAND SAW TO CUT THROUGH ONE ROW OF CUPS AS SHOWN. ALWAYS CUT PARALLEL TO RIBS, IN ORDER TO MAINTAIN MORE STRENGTH IN THE RACK.



0248 TRIMMING INSTRUCTIONS

***** NOTE *****

SOME INCUBATOR TRAYS ARE TOO SHORT TO ACCEPT TWO RACKS.

WHEN THIS HAPPENS, USE A KNIFE, HACK SAW OR BAND SAW TO CUT THROUGH ONE ROW OF CUPS AS SHOWN. ALWAYS CUT PARALLEL TO RIBS, IN ORDER TO MAINTAIN MORE STRENGTH IN THE RACK.

